

PRODUCT DATASHEET

# KG1400 Auto- Locking Lockset with Secondary Override

The KG1400 Auto-Locking lockset is the product for providing a dual means of secure entry to a patient or staff area.



## Terms and conditions:

To see our full terms and conditions, please visit:  
[kingswaygroupglobal.com/terms](https://kingswaygroupglobal.com/terms)

# Your Partner in Patient Safety

## Lockset Information

### Patient Function

Patients can access from the outside using the anti-ligature pull handle and key cylinder if issued a key. The door will automatically lock shut when closed. Patients can always exit from the inside by turning the anti-ligature thumb turn pull.

### Staff Function

Staff can always unlock from the outside using the key even if the thumb turn is held from the inside. The door will automatically lock when closed.

If the patient has jammed the keyway, the lock can be opened using the Kingsway Override Key (8mm Square or Oval secure).

### Components

The KG1400 Auto-Locking Lockset with Secondary Override comprises of:

- 1 x KG40 Ergo Grip Pull Handle
- 1 x KG84 Thumb Turn Pull
- 1 x KG131 Auto-Locking Lockcase
- 1 x KG176 Euro Escutcheon

### Handle Options

The KG40 Ergo Grip Pull Handle can be replaced with the following handles as per project requirement:

- KG30 Easy Grip Pull Handle
- KG60 Classic Grip Pull Handle
- KG80 Recessed Pull Handle.

## Guarantee

### 5 Year Guarantee\*

\*For faulty manufacture and not for damage

### Specification

KG1400 Auto-Locking Lockset With  
Secondary Override