

PRODUCT DATASHEET

KG1100 Bedroom Lockset

The KG1100 Bedroom lockset is the ideal solution to allow flexible access to a patient bedroom or other patient accessible areas.



Terms and conditions:

To see our full terms and conditions, please visit:
kingswaygroupglobal.com/en-au/terms

Your Partner in Patient Safety

Lockset Information

Patient Function

The patient can access from the outside using the anti-ligature pull handle and key cylinder (if issued a key) or the door will be open on the roller bolt. The patient can always exit from the inside by unlocking the deadbolt with the anti-ligature thumb turn pull.

Staff Function

Staff can always unlock the door from the outside using the key cylinder even if the thumb turn is held from the inside. Staff can also lock the door from the outside using the key.

Note: the KG1100 should not be specified on an outward opening door in a patient area as it does not have a secondary means of override - specify KG1200 Bedroom Lockset with Secondary Override.

Components

The KG1100 Bedroom Lockset comprises of:
1 x KG40 Ergo Grip Pull Handle
1 x KG82 Thumb Turn Pull
1 x KG130 Roller Bolt Deadlock Lockcase
1 x KG176 Euro Escutcheon

Handle Options

The KG40 Ergo Grip Pull Handle can be replaced with the following handles as per project requirement:

KG30 Easy Grip Pull Handle
KG60 Classic Grip Pull Handle
KG80 Recessed Pull Handle.

Warranty

1 Year Manufacturer's Defect Warranty*

*For faulty manufacture and not for damage

Specification

KG1100 Bedroom Lockset