

PRODUCT DATASHEET

# KG1200 Bedroom Lockset with Secondary Override

The KG1200 Bedroom lockset is the ideal solution to provide dual means of access to a patient bedroom or other patient accessible area.



## Terms and conditions:

To see our full terms and conditions, please visit:  
[kingswaygroupglobal.com/en-au/terms](https://kingswaygroupglobal.com/en-au/terms)

# Your Partner in Patient Safety

## Lockset Information

### Patient Function

The patient can access from the outside using the anti-ligature pull handle and key cylinder (if issued a key) or the door will be open on the roller bolt. The patient can always exit from the inside by unlocking the deadbolt with the anti-ligature thumb turn pull.

### Staff Function

Staff can unlock the door from the outside using the key cylinder even if the thumb turn is held from the inside. Staff can also lock the door from the outside using the key. If the patient has jammed the keyway, the deadbolt can be released by using the Kingsway Override Key (8mm Square or Oval secure).

### Components

The KG1200 Bedroom Lockset comprises of:  
1 x KG40 Ergo Grip Pull Handle  
1 x KG84 Thumb Turn Pull with External Override  
1 x KG130 Roller Bolt Deadlock Lockcase  
1 x KG176 Euro Escutcheon

### Handle Options

The KG40 Ergo Grip Pull Handle can be replaced with the following handles as per project requirement:

KG30 Easy Grip Pull Handle  
KG60 Classic Grip Pull Handle  
KG80 Recessed Pull Handle.

## Warranty

### 1 Year Manufacturer's Defect Warranty\*

\*For faulty manufacture and not for damage

### Specification

KG1200 Bedroom Lockset with Secondary Override